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Potentiometric Titrations are useful for determining a sample's IEP. As mentioned above, the IEP may be desirable or otherwise. Potentiometric Titration plots may display plateau regions for Zeta potential vs. pH. Such measurements enable manufacturers to optimize use of acids or bases for transportation and storage. Figure 2 shows a potentiometric titration on a Ludox-TM silica sample by automatic addition of 1N HCl. This titration was performed automatically by a Zeta-APS, Zeta Acoustic Particle Sizer, instrument from Matec Applied Sciences, Northborough, MA USA 3. Figure 2 shows that below pH 4, the Zeta curve approaches a plateau region while Conductivity increases more rapidly. This suggests that the silica particles are becoming saturated with H⁺ ions as the pH is lowered. Conductivity increases more rapidly as more of these H⁺ ions stay in the continuous water solvent as opposed to diffusing through the slipping plane toward the particle surface.

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
Figure 2. Automatic potentiometric titration of 10%-vol Ludox-TM by 1N HCl addition.

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